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N.A.

### 1. Goal

To describe how Peripheral Mononuclear Blood Cells (PBMC's) must be isolated and stored within the framework of the Parelshoer Institute for the purposes of the Radboud Biobank, so that:

1. Potential users know how the biomaterial has been handled.
2. The laboratory preparing the biomaterial for storage knows the minimal requirements for the handling of it.

#### 1.1 Scope off application

This procedure is applicable to all employees concerned at the Radboud Biobank and the Radboud Laboratory of Medical Immunology.

## 2. Protocol PBMC's Standard

### Collection

Delivery	in 2-4 collection tubes.
Type of tube	10 ml lithium-heparin tube. Standard tube, no gel, no protease inhibitors. BD tube 367526 or equivalent.
Temperature	Keep the tubes at room temperature until preparation.

### Preparation

Time until freezing	Within 24 hours after collection.
Isolation protocol	Isolation of cells from the blood by density gradient centrifugation.
Isolation medium	Density of separation medium must be 1,077 g/ml (Lymphoprep™ 1.077).

### Storage

Register	Register the microtubes in the Biobank management system and attach a cryovial label.
Storage medium	20% freezing medium: 6 parts RPMI, 2 parts inactivated flebopool and 2 parts DMSO. Add in 1:1- ratio to culture medium.
Aliquotation	Divide material over 2 microtubes of $5 \times 10^6$ – $20 \times 10^6$ cells. When there are more than $40 \times 10^6$ cells, divide over more microtubes. Register the exact number in the Biobank management system.
Type of tube	2 ml cryotube™ (Nunc)
Temperature	Freeze material first at $-80^{\circ}\text{C}$ (min. 24 hours, max. 1 week). Long-term storage in liquid nitrogen ( $-196^{\circ}\text{C}$ ).

### Miscellaneous

Deviations	Note all deviations from this procedure in the Biobank management system.
------------	---

### 3. Quality assurance

Review takes place via an internal audit (see the audit plan).

#### 3.1 Performance indicator

95% of the relevant samples have been prepared and stored within the parameters of this procedure.

### 4. Comments

This SOP fits within the framework of the Parelinoer Institute. Because the PSI standards are not unambiguous, some personal choices have been made concerning aliquotation. For explanation see document Bio-17165 "Onderbouwing SOP Isolatie en opslag van PBMC's".

### 5. Relevant documents

Qdocs, literature, legislation and codes of conduct, website etc.

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