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Revision changes

21-03-2019 Time until preparation is changed in time until freezing.

21-03-2019 REMP-tube is changed in MATRIX-tube

1. Goal

To describe how DNA from blood (500 µl) must be isolated and stored for the purposes of the Radboud Biobank, so that:

1. Potential users know how the biomaterial has been handled.
2. The laboratory preparing the biomaterial for storage knows the minimal requirements for the handling of it.

1.1 Scope of application

This procedure is applicable to all employees concerned at the Radboud Biobank and the Department of Human Genetics.

2. Protocol isolation and storage of DNA from blood (500 µl)

Collection

Delivery	Minimum of 500 µl EDTA whole blood
Type of needle	No requirements, is not considered critical.
Temperature	Keep at room temperature (RT) or at -80°C depending on the expected time until isolation takes place.

Preparation

Time until freezing	Isolate DNA within 48 hours of collection at RT or within a maximum of 4 months if stored at -80°C.
Isolation protocol	Isolate from the 500 µl EDTA blood received.
DNA-concentration and volume	After isolation determine the DNA-volume, DNA-concentration (in µg/ml) and the OD-ratio 260-280 nm from the stock solution and register in the Biobank management system.

Storage

Aliquotation	Store the biomaterial in a stock solution as 1 aliquot. N.B. The first time the DNA is issued make a working solution of 200 µl and 100 µg/ml and store in MATRIX-tube. After 5 issuances, make a new working solution from the stock solution at the time that the sixth issuance is requested.
Type of tube	MATRIX-tube
Temperature	-20°C

Miscellaneous

Data management	Register the following in the biobank management system: 1. Date blood was collected 2. Date and time of storage at -20°C (is date of DNA-isolation). 3. (pseudonymised) identification number, name of sub-biobank/ pearl/ project. 4. Volume, OD-ratio en DNA-concentration. 5. Deviations from the specified time until preparation. 6. Date and time of working solution is stored at -20°C (is the date the working solution is made), volume and concentration of the working solution. 6. Location of the DNA in the storage facility.
-----------------	---

3. Quality assurance

Review takes place via an internal audit (see the audit plan).

3.1 Performance indicator

95% of the sub biobank samples have been prepared and stored within the parameters of this procedure.

4. Comments

-

5. Relevant documents

Qdocs, literature, legislation and codes of conduct, website etc.

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