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Revision changes

06052019 Time until preparation is changed in time until freezing.
06052019 REMP-tube is change in MATRIX tube.

1. Goal

To describe how DNA from saliva must be isolated and stored for the purposes of the Radboud Biobank, so that:

1. Potential users know how the biomaterial has been handled.
2. The laboratory preparing the biomaterial for storage knows the minimal requirements for the handling of it.

1.1 Scope of application

This procedure is applicable to all employees concerned at the Radboud Biobank and the Department of Human Genetics.

2. Protocol isolation and storage of DNA from saliva

Collection

Preanalysis	Conditions for the collection of saliva for DNA investigations (Bio-17145).
Collection	2 ml saliva (or max. 5 “sponges” from Oragene) added to Oragene/saliva vial (Oragene kit; Oragene.DNA (OG-500) from DNAgenotek, Ontario, Canada).
Temperature	Keep tubes at room temperature until preparation.

Preparation

Time until freezing	One year maximum (shelf live of the saliva sample after combining with the Oragene-solution).
Isolation protocol	Isolate the DNA from the saliva (the whole Oragene/saliva vial) in accordance with the protocol “chemagic DNA Saliva Kit special” from DNA-genotek. Estimated yield of approximately 80 µg DNA.
DNA-concentration and volume	After isolation determine the DNA-volume, DNA-concentratie (in µg/ml) en OD-ratio 260-280 nm of the stock solution and register in the Biobank management system.

Storage

Aliquotation	Store the material in a stock solution in 1 aliquot. N.B. The first time the DNA is issued make a working solution of 200 µl and 100 µg/ml and store in MATRIX tube. After 5 issuances, make a new working solution from the stock solution at the time that the sixth issuance is requested.
Type of tube	MATRIX-tube
Temperature	-20°C

Miscellaneous

Data management	Register the following in the Biobank management system: 1. Date saliva was collected 2. Date and time of storage at -20°C (is date of DNA-isolation). 3. (pseudonymised) identification number, name of sub-biobank/ pearl/ project. 4. Volume, OD-ratio en DNA-concentration. 5. Deviations from the specified time until preparation. 6. Date and time working solution is stored at -20°C (is
-----------------	---

the date the working solution is made), volume and concentration of the working solution. 6. Location of the DNA in the storage facility.

3. Quality assurance

Review takes place via an internal audit (see the audit plan).

3.1 Performance indicator

95% of the sub biobank samples have been prepared and stored within the parameters of this procedure.

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