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## Revision changes

Comment added about the delivery of samples from the Biobank Ophthalmology and Pearl Hemophilia.

08042019 time until preparation is changed in time until freezing

08042019 REMP-tube is changed to MATRIX-tube

26-10-2020: Format is changed

26-10-2020: Removed: For Parelsnoer biobanks pass along the concentration an OD-ratio to the central database (Radboud Biobank employee task)

26-10-2020: Concentration in ng/ul instead of ug/ml

26-10-2020: Specifications Matrix tube added

26-10-2020: Removed: After 5 issuances, make a new working solution from the stock solution at the time that the sixth issuance is requested.

26-10-2020: Origin is added to data management

03-05-2021: making of working solution the first time the DNA is issued is removed from this SOP and is never done.

03-05-2021: Added: If DNA concentration less than 30 ng/ul the isolation is failed, see Deviations.

03-05-2021: Temperature is changed: Removed from SOP because it is not done: Store blood at -80°C when isolation is not possible within 48 hours. For the quality of the blood it is better to store longer at RT. Only batches blood are stored at -80 until isolation is planned.

03-05-2021: Time until freezing is changed

03-05-2021: Removed from because this is never checked: Isolate DNA from blood stored at -80 °C within a maximum of 4 months

## 1. Goal

To describe how DNA from blood must be prepared and stored for purposes of the Radboud Biobank (RB), so that:

1. Interested parties including potential users know how the biomaterial has been handled.
2. The laboratory can assess whether they can process the biomaterial as described under paragraph 2 "Preparation".
3. The laboratory knows the requirements for the collection, processing and storage of the biomaterial to enable registration of deviations.
4. The sub biobank or project knows the requirements for the collection and transport of the biomaterial to ensure correct delivery.
5. The RB is able to couple this procedure to the biomaterial in storage to determine the fitness for purpose.

### 1.1 Scope of application

This procedure is applicable to all employees concerned at the RB, the Sub-biobank/Project and the Department of Human Genetics of the Radboudumc.

## 2. Protocol DNA isolation from blood

### Collection

Delivery	In collection tube.
Type of tube	Standard EDTA tube, no gel and protease inhibitors. At least 4 ml; 6ml is preferred.
Type of needle	No requirements is not considered critical.
Temperature	Keep blood at room temperature (RT) until isolation. Only batches are stored at -80 °C until the isolation is planned.

**Title: SOP Isolation and storage of DNA from blood****Preparation**

Time until freezing	Isolate the DNA as soon as possible (blood stored at RT). Aim: within 48 hours Maximum: within 1 week
Isolation protocol	DNA isolation from EDTA-whole blood with ChemagicStar ( <a href="#">Qdoc 084294</a> ).
DNA-concentration and volume	After isolation determine the DNA-concentration (in ng/ul) and the OD-ratio 260-280 nm from the stock solution and register in the Biobank management system.

**Storage**

Aliquotation	Store the DNA in a 500 ul stock solution (1 aliquot).
Type of tube	MATRIX-tube (Matrix™ 2D Barcoded Open-Top Storage Tubes, ThermoFisher Scientific).
Temperature	-20°C (Liconic)

**Miscellaneous**

Data management	Register the following data in the biobank management system: <ol style="list-style-type: none"> <li>1. Date and time blood was collected</li> <li>2. Input-date is time of storage at -20°C after DNA-isolation.</li> <li>3. Identification number (pseudomised)</li> <li>4. Name of sub-biobank/ pearl/ project.</li> <li>5. Fraction number</li> <li>6. OD-ratio and DNA-concentration.</li> <li>7. Origin of the sample</li> </ol>
Deviations	Note all deviations from this procedure in the Biobank management system. If the DNA concentration is less than 30 ng/ul isolation has failed. Failed isolation is registered, and the RB is informed. The stocks are stored in the Liconic.

**3. Quality assurances**

Review takes place via i) an internal audit (see the [audit year plan](#) in DMS) and ii) 2-monthly performance reviews of deviations registered in the biobank management system (DMS map "[Registraties en controles](#)").

**3.1 Performance indicator**

95% of the sub-biobank samples have been prepared and stored within the parameters of this procedure.

**4. Comments**

1. The Ophthalmology Biobank delivers the remainder of a 10 ml EDTA-tube from which plasma has already been isolated. The Hemophilia Biobank uses the remainder of a 4- or 9-ml citrate-tube (instead of 6ml EDTA-tube).
2. Although this SOP is used for collection of Parelsnoer DNA, this SOP differs from the standard of Parelsnoer. Parelsnoer wants a minimum of 2 and prefers 3 or 4 aliquots. Given the costs and the frequency DNA is issued the RB decided to store only 1 DNA stock per isolation.

## 5. Relevant documents

Qdocs, literature, legislation and codes of conduct, website etc.

link	Title
<a href="#">Qdoc 083613</a>	Monstermanagement SOP KGCN
<a href="#">Qdoc 045505</a>	PVS 10.05 Primair laboratoriumproces, DNA-isolatiefaciliteit
<a href="#">Qdoc 072370</a>	DNA_NORM_ALLIN
<a href="#">Qdoc 084294</a>	DNA_isolatie_ALLIN
<a href="#">Qdoc 070512</a>	Liconic Bio.Li.X
<a href="#">Qdoc 083792</a>	Liconic vriesrobot, werking koelsysteem en alarmering
<a href="#">Qdoc 069283</a>	Bewaartermijnen Genetica
<a href="#">Qdoc 015298</a>	Bloedafname RLD: veneus en capillair

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**Title: SOP Isolation and storage of DNA from blood**

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