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Revision changes

01-04-2019 REMP tube is changed in MATRIX tube
 26-10-2020: Format is changed
 26-10-2020: Concentration in ng/ul instead of ug/ml
 26-10-2020: Specifications Matrix tube added
 26-10-2020: Removed: After 5 issuances, make a new work solution from the stock solution at the time that the sixth issuance is requested.
 26-10-2020: Origin is added to data management

1. Goal

To describe how DNA from citrate pellet must be prepared and stored for purposes of the Radboud Biobank, so that:

1. Interested parties including potential users know how the biomaterial has been handled.
2. The laboratory can assess whether they can process the biomaterial as described under paragraph 2 "Preparation".
3. The laboratory knows the requirements for the collection, processing and storage of the biomaterial to enable registration of deviations.
4. The sub biobank or project knows the requirements for the collection and transport of the biomaterial to ensure correct delivery.
5. The Radboud Biobank is able to couple this procedure to the biomaterial in storage to determine the fitness for purpose.

1.1 Scope of application

This procedure is applicable to all employees concerned at the Radboud Biobank, the Sub biobank/Project and the Department of Human Genetics of the Radboudumc.

2. Protocol DNA isolation from blood

Collection

Delivery	Citrate pellet form one citrate tube.
Type of tube	Standard 4,5 or 9 ml vacutainer citrate tube (Greiner or BD).
Type of needle	No requirements, is not considered critical.
Temperature	Keep at room temperature (RT) until isolation or store at -80°C when isolation is not possible within 48 hours.

Preparation

Time until freezing	Isolate DNA within 48 hours of collection, if blood is stored at -80 °C isolate DNA within a maximum of 4 months.
Isolation protocol	Isolate DNA from citrate pellet from the delivered citrate tube with ChemagicStar (Qdoc 047239).
DNA-concentration and volume	After isolation determine the DNA-concentration (in ng/ µl) and the OD-ratio 260-280 nm from the stock solution and register in the Biobank management system.

Storage

Aliquotation	Store the DNA in a 500 µl stock solution (1 aliquot). N.B. The first time the DNA is issued make a working solution of 200 µl and 100 ng/µl and store in MATRIX-tube.
--------------	--

Type of tube	MATRIX-tube (Matrix™ 2D Barcoded Open-Top Storage Tubes, ThermoFisher Scientific).
Temperature	-20°C

Miscellaneous

Data management	Register the following data in the biobank management system: 1. Date and time blood was collected 2. Input-date is time of storage at -20°C after DNA-isolation. 3. Identification number (pseudomised) 4. Name of sub-biobank/ pearl/ project. 5. Fraction number 6. OD-ratio and DNA-concentration. 7. Date and time of working solution is stored at -20°C (is the date the solution is made), volume and concentration of the working solution. 8. Location of the DNA in the storage facility. 9. Origin of DNA sample (Blood or Saliva)
Deviations	Note all deviations from this procedure in the Biobank management system

3. Quality assurances

Review takes place via i) an internal audit (see the [audit year plan](#) in DMS) and ii) 2-monthly performance reviews of deviations registered in the biobank management system (DMS map ["Registraties en controles"](#)).

3.1 Performance indicator

95% of the sub biobank samples have been prepared and stored within the parameters of this procedure.

4. Comments

5. Relevant documents

Qdocs, literature, legislation and codes of conduct, website etc.

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Title: SOP Isolation and storage of DNA from Citrate pellet**1. Doel**

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