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Revision changes

1. Goal

To describe how HPV smearstest material must be prepared and stored in the Radboud Biobank, so that:

1. Potential users know how the biomaterial has been handled.
2. The laboratory preparing the biomaterial for storage knows the minimal requirements for the handling of it.

1.1 Scope of application

This procedure is applicable to all employees concerned at the Radboud Biobank and the Radboud Pathology department.

2. Protocol HPV Smearstest

Collection

Delivery Batchwise each 3 months.

Type of tube 2 ml polypropylene microtube and screw-top with O-ring.
GREI722301UMC and GREI368380UMC

Temperature Keep tubes at room temperature.

Preparation

Preparation protocol Microtubes must be registered in the Biobank Management system.

Storage

Register Register the microtubes in the Biobank Management system.

Type of tube 2 ml polypropylene microtube and screw-top with O-ring.
GREI722301UMC and GREI368380UMC

Temperature Room temperature.

Comments The cell material is conserved in ThinPrep PreservCyt solution (Hologic Inc.). PreservCyt consists of approximately 50% methanol. Due to fire hazard and the toxicity of this storage solution these samples are stored in room M851. -1.221 (Storage bunker for hazardous substances, Central warehouse).

Miscellaneous

Deviations Note all deviations from this procedure in the Biobank Management system.

3. Quality assurance

Review takes place via an internal audit (see the audit plan).

3.1 Performance indicator

95% of the relevant samples have been prepared and stored within the parameters of this procedure.

4. Comments

When the cell material is issued it is transferred to the department of Pathology. The samples return to the Biobank because they can be issued 4 times on average.

5. Relevant documents

Qdocs, literature, legislation and codes of conduct, website etc.

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