How to send samples?

A centrifuged aliquot of cerebrospinal fluid (minimal volume 1 ml) can be sent to our laboratory for analysis.

Note: The CSF should be preferably sent frozen on dry ice to our laboratory; only when it can be delivered to our laboratory within 24 hours after withdrawal dispatch at room temperature is possible.

Analysis frequency

Every four weeks.

Turn-over time

Maximally 4.5 weeks.

Reference values

(Valid from 08-09-2017.)

Age	Upper reference limit (p95) (μg/L)
< 15 years	0.36
15 - 50 years	0.56
> 50 years	0.92

Address for sending material:

Radboud University Medical Center Radboudumc Diagnostic Laboratory Route 815 Attn. TML, Neurochemistry, dr. MM Verbeek PO 9101 6500 HB NIJMEGEN The Netherlands

Contact information : phone: #-31-24-3614567 fax: #31-24-3668754 e-mail: marcel.verbeek@radboudumc.nl

<u>Analysis costs</u> Please contact us