

PR007774	NEUROFILAMENT LIGHT CHAIN (NFL) PROTEINS IN CEREBROSPINAL FLUID (CSF)	PAGE 1/1 VALID UNTIL : 01-01-2024
----------	--	--------------------------------------

How to send samples?**Important!!!!**

Cerebrospinal fluid should be **collected, stored and shipped** in tubes made of polypropylene.

A centrifuged aliquot of cerebrospinal fluid (minimal volume 1 ml.; desirable volume 2 ml.) can be sent to our laboratory for analysis.

Note: The CSF should be preferably sent frozen on dry ice to our laboratory; only when it can be delivered to our laboratory within 24 hours after withdrawal dispatch at room temperature is possible.

Address for sending material:

Radboud University Medical Center
Radboudumc Diagnostic Laboratory
Route 815
Attn. TML, Neurochemistry, dr. MM Verbeek
PO 9101
6500 HB NIJMEGEN
The Netherlands

Contact information :

phone: #-31-24-3614567

fax: #31-24-3668754

e-mail: marcel.verbeek@radboudumc.nl

Analysis frequency

Every 4 weeks

Turn-over time

Maximally 5 weeks

Reference values

Age	NFL (ng/L)
< 15 yrs	< 1000
15 – 50 yrs	< 2200
> 50 yrs	< 2700

Analysis costs

Please contact us